ACETYLATION

Substitution of the active hydrogen of phenols, alcohols, thiols or amines (primary and secondary) by an acetyl group (CH₃CO) (Figure 1) is known as acetylation. Acetyl chloride, acetic anhydride or acetic acid can be used as acetylating agent.

• Acetyl chloride:

- **1.** It reacts vigorously.
- **2.** Liberating HCl which converts half of the amine to its hydrochloride salt rendering it incapable of participating in the reaction.
- 3. Has a tendency to get hydrolysed by moisture.
- 4. The reaction is not reversible.
- Acetic acid:
- 1. Its reaction with the compound is reversible and slow (requires long heating).
- Acetic anhydride: is preferred for a laboratory synthesis
- 1. Acetic anhydride is easy to handle
- **2.** The reaction is not reversible.
- 3. The procedure gives a product of high purity and in good yield.




X= NH, NR, S or O

Figure 1: General structure of an acetyl derivative.


Organic compound	Acetylation product	Organic synthesis purpose	Organic Pharmaceutical chemistry purpose
R-OH	Esters	Protection	Prodrugs, soft drugs
Ar-OH	Esters	Protection	Prodrugs, soft drugs
R-SH	Thioesters	Protection	Prodrugs, soft drugs
Ar-SH	Thioesters	Protection	Prodrugs, soft drugs
R-NHR*	Amides	Protection (not common)	Drugs
* R= H or R groups			

DRUGS AS ACETYL DERIVATIVES




CHEMISTRY OF ACETYLATION

Compounds with active hydrogen will react with acetic anhydride to produce acetyl derivatives by nucleophilic substitution reaction as in scheme (1):



The same mechanism will be operated if **O** or **S** was used instead of **N**



PROCEDURE

- a. Weigh 1.4g of salicylic acid and transfer to clean, dry round bottomed flask.
- **b.** Add 5ml of acetic anhydride and 5 drops of concentrated H₂SO₄, use stopper and shake manually for seconds (clear solution may be obtained).
- **c.** Float the round bottomed flask in water bath [heated to 60-80°C] for 10-15minutes with stirring.
- d. Discontinue heating and immediately deliver 2ml of DW, why?
- e. When the flask is sufficiently cool, remove the stopper, and add 20ml of DW, why?
- **f.** Immerse the flask in ice bath with scratching by glass rod until crystallization is completed.
- g. Isolate the product (ASA), dry, and weigh.



QUESTIONS

- The equipments must be dry in this experiment, why?
- Heating is by water bath and not direct heating; why?
- Sulfuric acid has two possible roles in this experiment, what are they?
- Acetylation occurred on phenolic OH not carboxylic OH; why?
- What are the purposes of acetylation in organic synthesis and organic pharmaceutical chemistry?
- Acetylation and acylation; any different between them??
- *How to check by chemical tests the progress of the following pathway:*

